

Brain Imaging Package




Upgrade to

MrT
3300S




PARAMED
MEDICAL SYSTEMS

MrT 3300S



Brain Coil



3D viewer

- **Brain Coil** - Multi-array coil (2 channels) dedicated to Brain Imaging
- **Enlarged FOV** - Tool for geometrical distortion correction allows to scan with a larger Field of View
- **Pulse sequences package** for Brain Imaging - Includes 2D and 3D Angio MR, T1 and T2 weighted FLAIR
- **3D viewer** with Post-processing capabilities - including
 - MIP (Maximum Intensity Projection)
 - MPR (Multi Planar Reconstruction)
- Enhanced patient table functionalities - including memory function for **pre** and **post contrast** imaging
- Set of **positioning** pads and accessories ergonomically designed for Brain examination



PARAMed
MEDICAL SYSTEMS

Corporate Headquarters
Paramed srl

Corso F.M. Perrone, 73r - 16152 Genova - Italy
Tel +39 010 7404530 • Fax +39 010 7404530

USA branch office
Paramed Medical Systems Inc.

39 High Street, North Andover, MA 01845 - USA
Tel +1 (978) 327-5853 • Fax +1 (978) 975-9930

www.paramedmedicalsystems.com

info@paramed.it